

# Archived Document

This archived document is no longer being reviewed through the CLSI Consensus Document Development Process. However, this document is technically valid as of January 2017. Because of its value to the laboratory community, it is being retained in CLSI's library.



March 2001

---

## AUTO05-A

### Laboratory Automation: Electromechanical Interfaces; Approved Standard

This document provides standards for the development of an electromechanical interface between instruments and specimen processing and handling devices used in automated laboratory testing procedures.

---

A standard for global application developed through the Clinical and Laboratory Standards Institute consensus process.

# Clinical and Laboratory Standards Institute

*Setting the standard for quality in medical laboratory testing around the world.*

The Clinical and Laboratory Standards Institute (CLSI) is a not-for-profit membership organization that brings together the varied perspectives and expertise of the worldwide laboratory community for the advancement of a common cause: to foster excellence in laboratory medicine by developing and implementing medical laboratory standards and guidelines that help laboratories fulfill their responsibilities with efficiency, effectiveness, and global applicability.

## Consensus Process

Consensus—the substantial agreement by materially affected, competent, and interested parties—is core to the development of all CLSI documents. It does not always connote unanimous agreement, but does mean that the participants in the development of a consensus document have considered and resolved all relevant objections and accept the resulting agreement.

## Commenting on Documents

CLSI documents undergo periodic evaluation and modification to keep pace with advancements in technologies, procedures, methods, and protocols affecting the laboratory or health care.

CLSI's consensus process depends on experts who volunteer to serve as contributing authors and/or as participants in the reviewing and commenting process. At the end of each comment period, the committee that developed the document is obligated to review all comments, respond in writing to all substantive comments, and revise the draft document as appropriate.

Comments on published CLSI documents are equally essential, and may be submitted by anyone, at any time, on any document. All comments are managed according to the consensus process by a committee of experts.

## Appeals Process

When it is believed that an objection has not been adequately considered and responded to, the process for appeals, documented in the CLSI Standards Development Policies and Processes, is followed.

All comments and responses submitted on draft and published documents are retained on file at CLSI and are available upon request.

## Get Involved—Volunteer!

Do you use CLSI documents in your workplace? Do you see room for improvement? Would you like to get involved in the revision process? Or maybe you see a need to develop a new document for an emerging technology? CLSI wants to hear from you. We are always looking for volunteers. By donating your time and talents to improve the standards that affect your own work, you will play an active role in improving public health across the globe.

For additional information on committee participation or to submit comments, contact CLSI.

Clinical and Laboratory Standards Institute  
950 West Valley Road, Suite 2500  
Wayne, PA 19087 USA  
P: +1.610.688.0100  
F: +1.610.688.0700  
[www.clsi.org](http://www.clsi.org)  
[standard@clsi.org](mailto:standard@clsi.org)

ISBN 1-56238-432-5  
ISSN 0273-3099

AUTO05-A  
Vol. 21 No. 5  
Replaces AUTO5-P  
Vol. 19 No. 23

---

# Laboratory Automation: Electromechanical Interfaces; Approved Standard

Volume 21 Number 5

Richard A. McPherson, M.D., Chairholder  
Raymond E. Karcher, Ph.D., Vice-Chairholder  
Mary F. Burritt, Ph.D.  
James Callaghan  
Norberto DeOliveira  
Michael W. LaCourt  
Matthew L. Mittino  
Takahiro Nagatomi

## Abstract

CLSI document AUTO05-A—*Laboratory Automation: Electromechanical Interfaces; Approved Standard* defines a standard connection between instruments and specimen processing and handling devices (including automated centrifuges, automated aliquoters, specimen integrity devices, and automated storage and retrieval systems) and automation systems. The user of the technology can thus create an optimally functioning automated laboratory environment. The issues addressed include:

- optimal configuration of the transport system in relation to the instrument;
- relationship of the instruments and/or specimen processing and handling devices to different types of transportation referenced to a defined point of reference (POR) (e.g., automated guided vehicles or “AGVs,” conveyors, and other mechanisms);
- responsibility for positioning of the specimen container in relation to the transport system and instrument;
- communication between the instruments and specimen processing and handling devices and the electromechanical interface transportation mechanism; and
- safety issues related to operation.

CLSI. *Laboratory Automation: Electromechanical Interfaces; Approved Standard*. CLSI document AUTO05-A (ISBN 1-56238-432-5). Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19087, USA 2001.

The Clinical and Laboratory Standards Institute consensus process, which is the mechanism for moving a document through two or more levels of review by the health care community, is an ongoing process. Users should expect revised editions of any given document. Because rapid changes in technology may affect the procedures, methods, and protocols in a standard or guideline, users should replace outdated editions with the current editions of CLSI documents. Current editions are listed in the CLSI catalog and posted on our website at [www.clsi.org](http://www.clsi.org). If your organization is not a member and would like to become one, and to request a copy of the catalog, contact us at: Telephone: 610.688.0100; Fax: 610.688.0700; E-Mail: [customerservice@clsi.org](mailto:customerservice@clsi.org); Website: [www.clsi.org](http://www.clsi.org).



Copyright ©2001 Clinical and Laboratory Standards Institute. Except as stated below, any reproduction of content from a CLSI copyrighted standard, guideline, companion product, or other material requires express written consent from CLSI. All rights reserved. Interested parties may send permission requests to [permissions@clsi.org](mailto:permissions@clsi.org).

CLSI hereby grants permission to each individual member or purchaser to make a single reproduction of this publication for use in its laboratory procedure manual at a single site. To request permission to use this publication in any other manner, e-mail [permissions@clsi.org](mailto:permissions@clsi.org).

### **Suggested Citation**

CLSI. *Laboratory Automation: Electromechanical Interfaces; Approved Standard*. CLSI document AUTO05-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2001.

### **Previous Edition:**

October 1999

### **Reaffirmed:**

March 2012

### **Archived:**

January 2017

ISBN 1-56238-432-5  
ISSN 0273-3099

## Contents

Abstract.....	v
Committee Membership.....	iii
Active Membership.....	v
Matrix of NCCLS Laboratory Automation Standards.....	xii
Preface to Laboratory Automation Standards.....	xv
Foreword.....	xix
1 Introduction.....	1
2 Scope.....	1
3 Definitions .....	1
4 System Functions.....	11
4.1 Laboratory Automation System (LAS) Functions.....	11
4.2 Instrument or Specimen Processing/Handling Device Functions.....	12
4.3 Instrument or Specimen Processing/Handling Device Vendor Responsibilities .....	12
4.4 Specimen Identification .....	13
4.5 Positioning and Adjustment Responsibilities .....	13
5 Physical Interface Specifications .....	14
5.1 Point of Reference (POR) .....	14
5.2 Relationship Between Specimen Container Position and Aspiration Probe: Target and Alignment Specifications.....	15
5.3 LAS Specifications .....	18
5.4 Clearance Zones.....	19
6 Standard Precautions.....	20
6.1 The Specimen .....	21
6.2 The Operator.....	21
6.3 The Environment .....	22
7 Human Factors.....	22
7.1 Human Factor Aids.....	23
7.2 Ergonomic Factors.....	23
References.....	24
Appendix. Analysis of Critical Tasks .....	26
Summary of Comments and Subcommittee Responses.....	27
Related NCCLS Publications.....	32

## Foreword

In 1996, NCCLS formed the Area Committee on Automation to develop global standards for automated laboratory systems that would result in automated system components that could be integrated into a laboratory, independent of manufacturer. Additionally, the area committee would foster communication of related issues and developments worldwide. These standards are developed to be flexible, such that they accommodate the needs of current systems as much as possible, but are targeted primarily toward future generations of laboratory automation equipment. Five interrelated standards were developed simultaneously to ensure that all aspects of automation were addressed as well as to ensure interconnectivity of the resulting designs.

AUTO05-A, *Laboratory Automation: Electromechanical Interfaces* was developed because the electromechanical interface between the instruments and specimen processing and handling devices and transportation systems is a critical part of the cost-effective manufacture, implementation, and operation of laboratory automation technology. The development of a standard, compatible connection between the instruments and specimen processing and handling devices and the automation systems should enable the user of the technology to create an automated laboratory environment that will function optimally for the user's individual laboratory.

In order to develop a standard that is truly global, input from all involved parties including instrument and specimen processing and handling device manufacturers, automation systems manufacturers and designers, government agencies, and laboratorians is necessary and has been sought. Representatives from all these arenas have participated in and contributed to the development and review of this standard. These individuals constitute the audience for this document, and we are grateful for their input. This has been an open process, and all viewpoints were considered valid and important.

These specifications are also intended to complement the interrelated NCCLS standards developed by other automation subcommittees and support overall operational goals for future development in laboratory instrumentation and automation:

AUTO1—*Laboratory Automation: Specimen Container/Specimen Carrier;*

AUTO2—*Laboratory Automation: Bar Codes for Specimen Container Identification;*

AUTO3—*Laboratory Automation: Communications with Automated Clinical Laboratory Systems, Instruments, Devices, and Information Systems;* and

AUTO4—*Laboratory Automation: Systems Operational Requirements, Characteristics, and Information Elements.*

## Key Words

Clearance zones, instrument or specimen processing and handling device, laboratory automation system (LAS), point of reference (POR)

# Laboratory Automation: Electromechanical Interfaces; Approved Standard

## 1 Introduction

Because the electromechanical interface between laboratory instruments and specimen processing and handling devices and transportation systems is a critical part of the cost-effective manufacture, implementation, and operation of laboratory automation technology, the NCCLS Area Committee on Automation and the Subcommittee on Electromechanical Interfaces developed this standard. The area committee and subcommittee believe that a standard, compatible connection between the instruments and specimen processing and handling devices and the automation systems will enable the user of the technology to create an optimally functioning automated laboratory environment.

## 2 Scope

Issues to be addressed in this standard include the optimal orientation of the transportation system in relation to the instrument and specimen processing and handling devices; relationship of the electromechanical interface to different types of transportation (e.g., automated guided vehicles, conveyors, and other transportation mechanisms) referenced to a defined point of reference; responsibility for positioning of the specimen container/specimen carrier in relation to the transportation system and instrument and specimen processing and handling devices; communication between the instruments and specimen processing and handling devices and the electromechanical interface/transportation mechanism; and safety and ergonomic issues related to operation.

Several issues were considered important to the proper functioning of the electromechanical interface, but were better addressed by one of the other interrelated standards. These issues include the placement of the bar-code label on the specimen container (AUTO2); specimen containers which would be supported (AUTO1); and communication between the laboratory automation system and the instrument and specimen processing and handling devices (AUTO3).

This standard fits into the series of interrelated NCCLS automation standards AUTO1—*Laboratory Automation: Specimen Container/Specimen Carrier*; AUTO2—*Laboratory Automation: Bar Codes for Specimen Container Identification*; AUTO3—*Laboratory Automation: Communications with Automated Clinical Laboratory Systems, Instruments, Devices, and Information Systems*; and AUTO4—*Laboratory Automation: Systems Operational Requirements, Characteristics, and Information Elements*.

## 3 Definitions<sup>b</sup>

Some of the computer-, automation-, or robotics-related terms used in the five interrelated NCCLS automation documents can be found in ANSI X3.172<sup>1</sup>, ANSI X3.182-1990<sup>2</sup>, ASTM D966<sup>3</sup>, ASTM E1013<sup>4</sup>, ASTM F149<sup>5</sup>, ASTM F1156<sup>6</sup>, IEEE 100,<sup>7</sup> IEEE 610,<sup>8</sup> IEEE 1007<sup>9</sup>, and HL7 Version 2.4<sup>10,11</sup>:

**ACK**, *n* – 1) A data field name for a general acknowledgment message as specified in the HL7 protocol (*HL7 V2.4*<sup>10</sup>); 2) A communication control character transmitted by a receiver as an affirmative response to a sender (*ASTM*).

**ADT**, *n* – 1) An abbreviation for admission, discharge, or transfer; 2) A data field in a hospital information system denoting admission, discharge, or transfer.

<sup>b</sup> Some of these definitions are found in NCCLS document NRSL8—*Terminology and Definitions for Use in NCCLS Documents*. For complete definitions and detailed source information, please refer to the most current edition of that document.

**Aliquot, *n*** – A portion of a specimen placed in a separate container to facilitate concurrent testing or to hold in reserve for future use; **NOTES:** a) The portion of the specimen is typically removed from the original specimen after initial processing, such as centrifugation, to obtain serum or plasma samples, and is considered to be chemically identical to all other subdivisions of an original sample of serum, plasma, urine, cerebral spinal fluid (CSF), etc.; b) It may be necessary to identify the aliquot as an individual specimen distinct from the original specimen in a collection container labeled with a unique identifier that may be linked to or associated with the primary collection container.

**Analyzer, *n*** – An instrument and/or specimen processing and handling device that performs measurements on patient specimens of quantitative, clinically relevant analytes; **NOTE:** A portion of a patient's specimen is consumed in the analytic process.

**ANSI, *n*** – Acronym/Initialism for American National Standards Institute; **NOTE: In Automation,** the Microsoft Windows ANSI character set is composed of ISO 8859/x plus additional characters.

**ASTM, *n*** – The official name of the organization formerly known as the American Society for Testing and Materials; **NOTE:** ASTM has developed various high- and low-level communications protocols.

**Audit trail, *n*** – An electronic log of transactions, detailing all events which have occurred in the laboratory automation system, including date and time of these events, which operator was responsible or directs processes, and any additional details.

**Automated, *adj*** – A characterization applied when all analytical processes, including sample and reagent uptake, sample/reagent interaction, chemical/biological analysis, result calculation, and result readout are mechanized. **NOTE:** These are usually controlled by a set of stored, modifiable instructions.

**Automated instrument, *n*** – A laboratory instrument that may or may not be connected to a laboratory information system (LIS), hospital information system (HIS), and/or laboratory automation system (LAS), which performs measurements on a patient's sample; **NOTE:** These instruments may have specific hardware and/or software modifications that allow interface to a laboratory automation system.

**Automation system, *n*** – An automation system refers to a variety of possible systems that can include some of the following types: automated instruments, laboratory information systems (LIS), laboratory automation systems (LAS), hospital information systems (HIS), and front-end processing devices.

**Bar code, *n*** – **1)** An array of parallel rectangular bars and spaces that creates a symbology representing a number or alphanumeric identifier; **2)** An array of rectangular lines and spaces that are arranged in a predetermined pattern following unambiguous rules and representing data that are referred to as characters (*ASTM F1156*<sup>6</sup>); **3)** An identification code consisting of a pattern of vertical bars whose width and spacing identifies the item marked; **NOTE:** The code is meant to be read by an optical input device, such as a bar-code scanner. Applications include retail product pricing labels, identification of library documents, and railroad boxcar identification. (*IEEE 610.2*<sup>8</sup>)

**Bar length, *n*** – The length of the bars in the bar code.

**Bottom of cap, *n*** – The farthest point from the top of the container/test tube that the cap reaches; **NOTE:** This point may be inside the tube.

**Bottom of container//Bottom of tube, *n*** – The portion of the container/test tube farthest from the cap (see **Point of reference**).

**Bottom of tube, *n*** – See **Bottom of container**.



**Related NCCLS Publications\***

- AUTO1**      **Laboratory Automation: Specimen Container/Specimen Carrier.** This document contains standards for design and manufacture of specimen containers and specimen carriers used for collection and processing of specimens, such as blood and urine, for testing on laboratory automation systems.
- AUTO2**      **Laboratory Automation: Bar Codes for Specimen Container Identification.** This document provides specifications for use of linear bar codes on specimen container tubes in the clinical laboratory and for use on laboratory automation systems
- AUTO3**      **Laboratory Automation: Communications with Automated Clinical Laboratory Systems, Instruments, Devices, and Information Systems.** This document provides standards to facilitate accurate and timely electronic exchange of data and information between the automated laboratory elements.
- AUTO4**      **Laboratory Automation: Systems Operational Requirements, Characteristics, and Information Elements.** This document provides standards of interest to operators for display of system status information such as specimen location, reagent supply, and warnings and alerts to support laboratory automation operations.
- GP2-A2-C**    **NCCLS Procedure Manual Template.** This computer template enables laboratorians to prepare consistent technical procedures in the NCCLS format. The template and its user manual, used along with the GP2-A3 guideline, provide a procedure format that is as easy to use as a word processing program. Procedures can be stored as individual files for easy retrieval and updating, or they can be networked through the local computer system for electronic distribution throughout the laboratory. The template format consists of tables for recording essential information for all procedures and an outline of key headings for incorporating procedure-specific details.
- GP2-A3**      **Clinical Laboratory Technical Procedure Manuals – Third Edition; Approved Guideline (1996).** This document provides guidance for the patient-testing community by addressing the design, preparation, maintenance, and use of paper or electronic technical procedure manuals.
- GP18-A**      **Laboratory Design; Approved Guideline (1998).** This guideline provides a foundation of information about laboratory design elements that can be used to help define the issues being considered when designing a laboratory.
- GP19-A**      **Laboratory Instruments and Data Management Systems: Design of Software User Interfaces and End-User Software Systems Validation, Operation, and Monitoring; Approved Guideline (1995).** This document identifies important factors that designers and laboratory managers should consider when developing new software-driven systems and selecting software user interfaces. Also included are simple rules to help prepare validation protocols for assessing the functionality and dependability of software.

---

\* Proposed- and tentative-level documents are being advanced through the NCCLS consensus process; therefore, readers should refer to the most recent editions.

**Related NCCLS Publications (Continued)**

**H18-A2**      **Procedures for the Handling and Processing of Blood Specimens; Approved GuidelineXSecond Edition (1999).** This guideline addresses multiple factors associated with handling and processing of specimens, and factors that can introduce imprecision or systematic bias into results.

**H38-P**      **Calibration and Quality Control of Automated Hematology Analyzers; Proposed Standard (1999).** This document addresses calibration and quality control strategies for multichannel hematology analyzers; assignment of values to calibrator materials; calibration using stabilized blood controls; internal quality control; pair difference analysis; and use of the weighted moving average (-/+<sub>B</sub>) method.

**M29-A**      **Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue; Approved Guideline (1997).** A consolidation of M29-T2 and I17-P, this document provides guidance on the risk of transmission of hepatitis viruses and human immunodeficiency viruses in any laboratory setting; specific precautions for preventing the laboratory transmission of blood-borne infection from laboratory instruments and materials; and recommendations for the management of blood-borne exposure.

SAMPLE

# Explore the Latest Offerings From CLSI!

As we continue to set the global standard for quality in laboratory testing, we are adding products and programs to bring even more value to our members and customers.



By becoming a CLSI member, your laboratory will join 1,600+ other influential organizations all working together to further CLSI's efforts to improve health care outcomes. You can play an active role in raising global laboratory testing standards—in your laboratory, and around the world.

Find out which membership option is best for you at [www.clsi.org/membership](http://www.clsi.org/membership).



Find what your laboratory needs to succeed! CLSI U provides convenient, cost-effective continuing education and training resources to help you advance your professional development. We have a variety of easy-to-use, online educational resources that make eLearning stress-free and convenient for you and your staff.

See our current educational offerings at [www.clsi.org/education](http://www.clsi.org/education).



When laboratory testing quality is critical, standards are needed and there is no time to waste. eCLIPSE™ Ultimate Access, our cloud-based online portal of the complete library of CLSI standards, makes it easy to quickly find the CLSI resources you need.

Learn more and purchase eCLIPSE at [clsi.org/eCLIPSE](http://clsi.org/eCLIPSE).

For more information, visit [www.clsi.org](http://www.clsi.org) today.

SAMPLE



CLINICAL AND  
LABORATORY  
STANDARDS  
INSTITUTE®

---

950 West Valley Road, Suite 2500, Wayne, PA 19087 USA

ISBN 1-56238-432-5

P: +1.610.688.0100 Toll Free (US): 877.447.1888 F: +1.610.688.0700

E: [customerservice@clsi.org](mailto:customerservice@clsi.org) [www.clsi.org](http://www.clsi.org)